



Standard Guide for Evaluation of Cleanroom Disinfectants¹

This standard is issued under the fixed designation E 2614; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ε) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This guide identifies important factors to consider when selecting a disinfectant for use in a cleanroom or similar controlled environment and recommends test methods suitable for evaluating disinfectants. The proper selection of disinfecting agents combined with in-house qualification testing is a key element to a successful disinfection program. Recent publications of regulatory guidance/standards such as Chapter <1072>, “Disinfectants and Antiseptics” in the United States Pharmacopoeia address the issue of disinfectant testing but there is very little published guidance on what criteria and test methods should be used for selection and efficacy testing of disinfectants that will be used in cleanrooms and controlled environments.

1.2 An understanding of microbiology and microbiological techniques is essential. Knowledge in the following areas is recommended: microorganisms, antimicrobial products (disinfectants, sporicides, and decontamination agents), and the chemistry of disinfection, mechanism of activity of disinfectants on cells, application procedures, cleanroom surfaces, and environmental conditions within a cleanroom. This information is available in several published texts listed in the bibliography.

1.3 The theoretical basis for disinfectant activity is not addressed in this guide. An understanding of the effect of disinfectant concentration on microbial reduction (concentration exponent) and kinetics is desirable in determining the use-dilution of different disinfectants and in using dilution to neutralize a disinfectant for efficacy testing. USP chapter <1072> provides further information on this topic.

1.4 Although the information in this guide is written for the cleanroom environment, many of the principles outlined in this standard are applicable to manufacturing/processing environments outside of the cleanroom.

1.5 Evaluation of disinfectants for biofilm control is outside the scope of this document. The reader is referred to Guide E 1427.

1.6 The values state in inch-pound units are to be regarded as standard. No other units of measurement are included in this standard.

2. Referenced Documents

2.1 ASTM Standards:²

E 1427 Guide for Selecting Test Methods to Determine the Effectiveness of Antimicrobial Agents and Other Chemicals for the Prevention, Inactivation and Removal of Biofilm

E 2111 Quantitative Carrier Test Method to Evaluate the Bactericidal, Fungicidal, Mycobactericidal, and Sporocidal Potencies of Liquid Chemical Microbicides

E 2197 Quantitative Disk Carrier Test Method for Determining the Bactericidal, Virucidal, Fungicidal, Mycobactericidal and Sporocidal Activities of Liquid Chemical Germicides

E 2315 Guide for Assessment of Antimicrobial Activity Using a Time-Kill Procedure

2.2 Other Standard:

United States Pharmacopeia 30, Chapter <1072> Disinfectants and Antiseptics, May 1, 2007³

AOAC, Chapter 6, 15th Edition Official Methods of Analysis of AOAC International, Chapter 6: Disinfectants, 15th ed, 1990⁴

AOAC Chapter 6, 17th Edition Official Methods of Analysis of AOAC International, Chapter 6: Disinfectants, 17th ed., 2000⁴

BS EN 1040 Chemical Disinfectants and Antiseptics—Basic Bactericidal Activity⁵

BS EN 1276 Chemical disinfectants and antiseptics—Quantitative Suspension Test for the Evaluation of Bactericidal Activity of Chemical Disinfectants and Antiseptics

¹ This specification is under the jurisdiction of ASTM Committee E35 on Pesticides and Alternative Control Agents and is under the direct responsibility of Subcommittee E35.15 Antimicrobial Agents.

Current edition approved April 1, 2008. Published May 2008.

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

³ Available from U.S. Pharmacopeia (USP), 12601 Twinbrook Pkwy., Rockville, MD 20852-1790, <http://www.usp.org>.

⁴ Available from AOAC International, 481 North Frederick Ave., Suite 500, Gaithersburg, Maryland 20877-2417, <http://www.aoac.org>.

⁵ Available from British Standards Institute (BSI), 389 Chiswick High Rd., London W4 4AL, U.K., <http://www.bsi-global.com>.

used in Food, Industrial, Domestic, and Institutional Areas⁵

BS EN 13697 Chemical Disinfectants and Antiseptics—Quantitative Non-Porous Surface Test for the Evaluation of Bactericidal and/or Fungicidal Activity of Chemical Disinfectants used in Food, Industrial, Domestic and Institutional Areas⁵

3. Terminology

3.1 Definitions:

3.1.1 *antimicrobial, n*—describes an agent that kills bacteria or suppresses their growth or reproduction.

3.1.2 *bioburden, n*—the number and type of viable microorganisms that can be estimated using prescribed recovery procedures.

3.1.3 *biocide, n*—a physical or chemical agent that kills organisms.

3.1.4 *carrier, n*—a surrogate surface or matrix that facilitates the interaction of test microorganisms and treatment.

3.1.5 *cleanroom, n*—an area equipped with HEPA filtered air and used for aseptic processing.

3.1.6 *contact time, n*—predetermined time that a test microorganism is exposed to the activity of a test material.

3.1.7 *disinfectant, n*—a physical or chemical agent or process that destroys pathogenic or potentially pathogenic microorganisms on inanimate surfaces or objects.

3.1.8 *efficacy, n*—the proven performance of a product established under defined conditions.

3.1.9 *effectiveness, n*—a measure of the performance of a product.

3.1.10 *inoculum, n*—the viable microorganisms used to contaminate a sample, device, or surface, often expressed as to number and type.

3.1.11 *neutralization, n*—the process for inactivating or quenching the activity of a microbiocide, often achieved through physical (for example, filtration or dilution) or chemical means.

3.1.12 *qualification, n*—to determine effectiveness in the context of a given process.

3.1.13 *sanitizer, n*—chemical or physical agent(s) used to reduce the number of microorganisms to a level judged to be appropriate for a defined purpose and/or claim.

3.1.14 *soil load, n*—chemical or physical material(s) included in a test procedure to simulate conditions or use.

3.1.15 *sporicide, n*—chemical or physical agent(s) that kill spores.

3.1.16 *substrate, n*—surface on which an organism can grow

3.1.17 *surfactant, n*—synthetic detergent

3.2 Acronyms:

3.2.1 *MSDS*—Material Safety Data Sheet

4. Summary of Guide

4.1 Selecting and qualifying the appropriate disinfection agents is an integral factor in developing a compliant cleaning and disinfection program for a cleanroom. Significant factors to consider when selecting disinfectants for use in cleanrooms and controlled environments are discussed in this guide. A

summary of the most common test methods used to determine disinfectant effectiveness is also presented.

5. Significance and Use

5.1 Requirements for aseptic processing areas include readily cleanable floors, walls, and ceilings that have smooth, non-porous surfaces; particulate, temperature, and humidity controls; and cleaning and disinfecting procedures to produce and maintain aseptic conditions. These conditions combined with careful and thorough evaluation of the chemical agents used for the cleaning and disinfection program should lead to achieving the specified cleanliness standards and control of microbial contamination of products. In recent years, the use of disinfectants in pharmaceutical, biotechnology, medical device facilities, and associated controlled environments has been the subject of scrutiny by regulatory agencies.

5.2 An effective cleaning and disinfection program in aseptic processing areas of a Good Manufacturing Practice (GMP) facility is critical to assure the quality of the products. Manufacturers are being held to a high standard when it comes to product sterility and regulatory agencies are increasingly asking for validation data to support sanitization and disinfection procedures. Regulatory authorities now expect evidence of the efficacy of disinfection agents against environmental isolates. The FDA Guideline for Aseptic Processing states, “the suitability, efficacy, and limitations of disinfecting agents and procedures should be assessed. The effectiveness of these disinfectants and procedures should be measured by their ability to ensure that potential contaminants are adequately removed from surfaces.”⁶

5.3 Basic knowledge regarding the effectiveness of different chemical agents against vegetative bacteria, fungi, and spores will aid in selecting chemical agents.

5.4 A good understanding of test methods used to assess disinfectant effectiveness is important. Most methods are adaptable allowing the user to customize the methods to their specific requirements.

6. Procedure

6.1 Selection Procedure:

6.1.1 In pharmaceutical, biotechnology, and associated industries, the selection of one sanitizer, one or two disinfectants, and one sporicide is typical. More than one type of chemical agent is needed to obtain the proper balance of effective microbial control and minimal surface damage because products vary in spectrum of activity and formulation. Currently, most facilities select one or two disinfecting agents to use on a routine basis and supplement with a sporicide which is used on a less frequent routine basis to address spores that may not be destroyed by the routine disinfectant. All three product types should be evaluated for effectiveness with appropriate test methods. Typical agents used in cleanrooms are:

⁶ FDA Guideline for Industry: Sterile Drug Products Produced by Aseptic Processing—Current Good Manufacturing Process, September 2004. Available at <http://www.fda.gov/cder/guidance/index.htm>

Sanitizers:	70 % v/v Isopropyl Alcohol (IPA) 70 % v/v Ethanol
Disinfectants:	Phenols Quaternary Ammonium Compounds Hydrogen Peroxide 3 % Sodium Hypochlorite < 0.10 % 70 % v/v Isopropyl Alcohol (IPA) 70 % v/v Ethanol
Sporicides:	Sodium Hypochlorite > 0.3 % Hydrogen Peroxide > 6 % Peracetic Acid Chlorine Dioxide Vaporized Hydrogen Peroxide

6.1.2 Sanitizers play a role in a cleanroom disinfection program. Isopropyl alcohol has good broad-spectrum efficacy, but in the manner which it is used routinely, minimal product contact time and flammability are of concern. Classifying 70 % IPA as a decontamination or residue control agent is more appropriate since it is often used during cleaning operations, as a sanitizer during manufacturing, and for removal of residues left from other antimicrobial agents.

6.1.3 It is important to choose a disinfectant that incorporates a surfactant in its formulation to help clean surfaces. Occasionally, the choice will be made to rotate two different disinfectants. These agents are designed for frequent use.

6.1.4 Sporocidal formulations typically contain highly reactive chemicals (for example, sodium hypochlorite or oxidizers) and may cause damage to surfaces, even stainless steel, if used too frequently. However, they are a necessity in the cleanroom since routine disinfectants are not effective against bacterial endospores or resistant mold spores such as *Aspergillus niger*. Typically, they are not used on a daily basis to avoid damage to substrates and reduce the potential for irritation to cleaning personnel.

6.1.5 Critical Factors in Selection of Chemical Agents:

6.1.5.1 *Product Information*—One of the first steps in selection of a chemical agent is to collect all pertinent and available information from the potential supplier. This would include EPA registration, technical data sheets, MSDS, recommended directions for use as well as data on substrate compatibility, stability, and microbiocidal efficacy performed according to acceptable standards.

6.1.5.2 *Cleanroom Bioburden*—A review of environmental monitoring data is necessary to determine the number and types of microorganisms to be controlled. Not all chemical agents are effective against all microorganisms and not all are sporocidal. Knowing the types of microbial flora in the area to be disinfected aids in selecting an appropriate disinfectant with the desired spectrum of activity. From this information, a facility can determine what efficacy is needed and identify potential product candidates.

6.1.5.3 *Contact Time*—For a disinfectant to be effective, a surface must remain wet for a defined period of time. The contact time may vary depending on the target microbial contamination, air temperature, as well as the type of surface being treated.

6.1.5.4 *Preparation*—Products can be purchased in concentrate form or a ready-to-use form. If a concentrated product is selected, dilution will be necessary. The accuracy of prepara-

tion of a use-dilution is critical for optimal performance. Regulatory standards require that disinfectants used in cleanrooms are sterile, therefore it may be necessary to purchase the agent sterile or sterilize it in-house. The kind of water available for dilution is also important. Water hardness or softness can interfere with the action of some disinfectants. It is recommended that purified water or better be used to dilute the concentrate.

6.1.5.5 *Safety and Toxicity*—The product MSDS should be reviewed to determine special storage and handling requirements and to identify any potential hazards to personnel working with the agent. Appropriate regulations should be reviewed regarding disposal.

6.1.5.6 *Compatibility*—Some disinfectants are incompatible with others and can cause problematic residues. Performance of disinfectants can be decreased if they are incompatible with cleaning agents used in the area. Additionally, compatibility with surface types in the cleanroom should be evaluated.

7. Disinfectant Challenge Testing

7.1 Once the chemical agents to be used in the cleanroom have been selected, they must be qualified for use with efficacy testing. This testing provides a means of assessing that a disinfectant will render a surface microbiologically safe for the intended purpose. A combination of several test methods is necessary to qualify a disinfectant for use.

7.1.1 Data from *in vitro* testing contributes to the overall picture of how a disinfectant will potentially perform in practice. *In vitro* testing includes suspension tests and carrier tests.

7.1.2 *In situ* testing demonstrates that the disinfectant, as used, will kill the target organisms.

7.1.3 Trend analysis of environmental monitoring data provides an on-going analysis on the effectiveness of the disinfectant.

7.2 Laboratory tests for the evaluation of disinfectants must be scientifically based and ensure that the agent is safe and effective for its intended use.

7.3 Following the same basic methods when evaluating several potential disinfectants for use is a good idea as it allows direct comparison between products.

7.4 Regulatory guidance documents have not identified specific methods for disinfectant qualification testing. Test protocols can be based on ASTM standards, the Association of Official Analytical Chemists (AOAC) methods, or the European standards (EN) for disinfectant efficacy testing. However, some modification of these methods is necessary when qualifying disinfectants for cleanrooms. These modifications include the use of environmental organisms, contact time, clean/dirty conditions, carrier type/material and preparation of the disinfectant.

7.4.1 Use of environmental organisms. Reference organisms such as *Escherichia coli*, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Candida albicans*, and *Aspergillus niger* should be included in the testing for microbiocides to be used in cleanrooms. However, regulatory authorities have been clear and consistent regarding microorganisms used for disinfectant qualification and expect environmental isolates from the cleanroom environment to be included in the testing as well.

Regulatory guidelines state that, “Routinely used disinfectants should be effective against the normal microbial vegetative flora recovered from the facility, “which substantiates including environmental isolates in addition to reference strains.⁶ It is not necessary to include all environmental isolates in disinfectant efficacy testing but it is necessary to have a valid rationale for the selection of the organisms that are included in the studies. Since these are the microorganisms most likely to present contamination risks to the product, it is critical to prove that the disinfectants used in a microbial control program can control them.

7.4.2 Contact Time—The contact time used for the test should be realistic. Regulatory authorities have cited facilities that fail to qualify disinfectants in the manner in which they are used, including disinfectant exposure time.⁷

7.4.3 Clean/Dirty Conditions—Some of the published methods test the disinfectant under “clean” and “dirty” conditions. The dirty condition represents an area known to contain levels of organic and inorganic materials. Generally, testing under clean conditions only is sufficient as most sterile manufacturing areas would have little, if any, organic or inorganic contamination.

7.4.4 Preparation—The disinfectant used for testing must be prepared the same way it will be used in practice. While some of the published AOAC and EN methods use hard water, this would not be the case for disinfectants used in cleanrooms.

7.5 Test Methods:

7.5.1 Suspension Tests—These tests are used to screen disinfectants for their efficacy at various concentrations and contact times against a wide range of environmental isolates and reference strains. They can be used for preliminary screening of several disinfectant products to determine which products should be tested further. The procedure involves pipeting a known volume of an organism directly into a test tube containing the sanitizer, disinfectant, or sporicide. After the required contact time, serial dilutions are performed in neutralization medium for enumeration. Subsequent plating of the positive and negative controls together with the neutralization confirmation and the test sample are completed. The following published methods describe methodology for suspension testing:

7.5.1.1 Guide E 2315,

7.5.1.2 European Standards including BS EN 1040 and BS EN 1276, and

7.5.1.3 AOAC Use-Dilution Method.

7.5.2 Surface Challenge Tests—This guide provides information about the performance of the chemical agent under laboratory conditions representative of practical use. Test organisms are dried on coupons (typically 2 by 2 in.) made of varying material substrates from the cleanroom (stainless steel, aluminum, epoxy, vinyl, plastic, terrazzo, and others). Once dried, the entire coupon is submerged in a test tube containing the disinfectant and subsequently (after the specified contact time) a serial dilution is performed in the same manner as the suspension tests. Surface testing presents a more stringent challenge to disinfectant products and is a closer representation

of practical in-use conditions. USP <1072> states that sufficient organisms need to be inoculated on a 2 by 2 in. square of the surface being decontaminated (that is, coupon) to demonstrate at least a 2 log reduction for bacterial spores and 3 log reduction for vegetative bacteria. The following published methods describe methodology for surface challenge testing:

7.5.2.1 Quantitative Carrier Test Method E 2111,

7.5.2.2 Quantitative Disk Test Method E 2197,

7.5.2.3 BS EN 13697, and

7.5.2.4 AOAC Hard Surface Carrier Test Method

7.5.3 In Situ Tests—In situ testing is an important phase of the qualification process but not easily carried out since it involves shutdown of the cleanroom. In situ evaluations are intended to provide data on the performance of the disinfectant under actual use conditions and in a way that clearly demonstrates performance. Regulatory authorities have cited firms for failure to perform evaluation under actual use conditions.⁷ Conducting the in situ evaluation under worst-case conditions is the ideal situation. As an example, in situ qualification can be scheduled around a preventative maintenance shutdown of the cleanroom. Sampling can be done during the shutdown (prior to disinfectant application) and again after the cleanroom is back in operation and disinfected according to procedure. Data review is conducted to determine if the disinfectant effectively reduced contamination in the sampled areas. A statistical comparison of the frequency of isolation and numbers of microorganisms isolated prior to and after the implementation of the new disinfectant is performed.

7.5.4 Trend Analysis of Environmental Monitoring Data—Environmental monitoring data should be trended to determine whether or not a disinfectant is effective. Frequent review of environmental monitoring data is necessary in order to assess the effectiveness of a selected disinfectant against the cleanroom bioburden. Environmental monitoring data over an entire year should be used for the comparison to allow for seasonal fluctuations.

7.6 Requalification:

7.6.1 It is suggested that requalification of a disinfectant would be required only if one or a combination of the following were to occur:

7.6.1.1 A change in the disinfectant formulation

7.6.1.2 A significant change in the number and/or types of organisms isolated from the environment

7.6.1.3 A modification of the surfaces in the cleanroom

8. Summary

8.1 Disinfectant qualification is expensive and time consuming but expected by regulatory authorities. However, the choice of disinfectants for contamination control is an important part of assuring the effectiveness of a cleaning and disinfection program for a cleanroom. Many factors must be considered when selecting and testing disinfectants for use in a cleanroom. Companies should consider all the relevant data and make a reasonable justification to support their individual disinfection program and product choice.

9. Keywords

9.1 carrier tests; cleanroom; disinfectant; disinfectant qualification; efficacy; in situ testing; suspension tests

⁷ GMP Trends, Issue No. 666, October 15, 2004.

BIBLIOGRAPHY

- (1) Block, S. S. ed., *Disinfection, Sterilization, and Preservation*, Lippincott, Williams and Wilkins, Philadelphia, PA, 2001.
- (2) Denyer, S., Hodges, N., and Gorman, S. eds., *Hugo and Russell's Pharmaceutical Microbiology*, Blackwell Science Ltd, Oxford, UK, 2004.
- (3) Jiminez, L., ed., *Microbial Contamination Control in the Pharmaceutical Industry*, Marcel Dekker, Inc., New York, NY., 2004.
- (4) Martinez, J. E. "What is disinfectant validation?," *Pharmaceutical Technology*, March 2006, pp.72–80.
- (5) Moldenhauer, J. ed., *Laboratory Validation—A Practitioner's Guide*. PDA, Bethesda, MD, 2003.
- (6) Sartain, E. "The ABCs of Disinfectant Validation," *CleanRooms*, March 2005.
- (7) Scott, D. and Alexander, B. "A Practical Approach to Disinfectant Tests," *Cleanroom Technology*, May 2007.

ASTM International takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.

This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, at the address shown below.

This standard is copyrighted by ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, United States. Individual reprints (single or multiple copies) of this standard may be obtained by contacting ASTM at the above address or at 610-832-9585 (phone), 610-832-9555 (fax), or service@astm.org (e-mail); or through the ASTM website (www.astm.org).