



Standard Test Methods for Determination of Manganese (IV) in Manganese Ores by Redox Titrimetry¹

This standard is issued under the fixed designation E465; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ε) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 These test methods cover the determination of manganese dioxide in amounts commonly found in manganese ore. The determination measures the amount of manganese (IV) present in the sample. The result may be expressed as available oxygen or as manganese dioxide. The following test methods are included and may be used interchangeably:

Sections

Test Method A (Ferrous Ammonium Sulfate)	9-13
Test Method B (Periodate (Sodium Oxalate))	14-18

1.2 The values stated in SI units are to be regarded as standard. No other units of measurement are included in this standard.

1.3 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.*

2. Referenced Documents

2.1 ASTM Standards:²

D1193 Specification for Reagent Water

E50 Practices for Apparatus, Reagents, and Safety Considerations for Chemical Analysis of Metals, Ores, and Related Materials

E135 Terminology Relating to Analytical Chemistry for Metals, Ores, and Related Materials

E882 Guide for Accountability and Quality Control in the Chemical Analysis Laboratory

¹ These test methods are under the jurisdiction of ASTM Committee **E01** on Analytical Chemistry for Metals, Ores, and Related Materials and are the direct responsibility of Subcommittee **E01.02** on Ores, Concentrates, and Related Metallurgical Materials.

Current edition approved Sept. 1, 2011. Published September 2011. Originally approved in 1972. Last previous edition approved in 2005 as E465 – 00 (2005). DOI: 10.1520/E0465-11.

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

3. Terminology

3.1 *Definitions*—For definitions of terms used in this test method, refer to Terminology **E135**.

4. Significance and Use

4.1 This test method is intended to be used for compliance with compositional specifications for manganese dioxide content in manganese ores. It is assumed that all who use these procedures will be trained analysts capable of performing common laboratory procedures skillfully and safely. It is expected that work will be performed in a properly equipped laboratory and that proper waste disposal procedures will be followed. Appropriate quality control practices shall be followed, such as those described in Guide **E882**.

5. Interferences

5.1 The elements ordinarily present in manganese ores do not interfere in either test method.

6. Reagents and Materials

6.1 *Purity of Reagents*—Reagent grade chemicals shall be used in all tests. Unless otherwise indicated, it is intended that all reagents conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society where such specifications are available.³ Other grades may be used, provided it is first ascertained that the reagent is of sufficient high purity to permit its use without lessening the accuracy of the determination.

6.2 *Purity of Water*—Unless otherwise indicated, references to water shall be understood to mean reagent water conforming to Type I or II of Specification **D1193**. Type III or IV may be used if they effect no measurable change in the blank or sample.

³ *Reagent Chemical, American Chemical Society Specifications*, American Chemical Society, Washington, DC. For suggestions on the testing of reagents not listed by the American Chemical Society, see *Analar Standards for Laboratory Chemicals*, BDH Ltd., Poole Dorset, U. K., and the *United States Pharmacopeia and National Formulary*, U.S. Pharmaceutical Convention, Inc., (USPC), Rockville, MD.

7. Hazards

7.1 For precautions to be observed in these methods, refer to Practice E50.

8. Sampling and Sample Preparation

8.1 The sample shall pass a No. 100 (150-µm) sieve.

TEST METHOD A—FERROUS AMMONIUM SULFATE METHOD

9. Summary of Test Method

9.1 The test sample is dissolved in an excess of ferrous ammonium sulfate solution. The manganese dioxide reacts with an equivalent amount of ferrous iron. The excess ferrous iron is titrated with standard potassium dichromate solution using sodium diphenylamine sulfonate as an indicator.

10. Reagents and Materials

10.1 *Ferrous Ammonium Sulfate Solution* (45 g/L)—Dissolve 45 g of ferrous ammonium sulfate [$\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$] in 1 L of sulfuric acid (H_2SO_4 , 1 + 7).

10.2 *Potassium Dichromate, Standard Solution* (0.1 N)

10.2.1 Dissolve 4.9035 g of Primary Standard Grade potassium dichromate ($\text{K}_2\text{Cr}_2\text{O}_7$) in water, transfer to a 1-L volumetric flask, dilute to volume, and mix.

10.3 *Sodium Diphenylamine Sulfonate Indicator Solution*, (2 g/L).

10.3.1 Dissolve 0.20 g of sodium diphenylamine sulfonate in 100 mL of water. Store in a dark-colored bottle.

11. Procedure

11.1 Transfer approximately 0.4 g of test sample to a small, dry weighing bottle and place into a drying oven. Dry at 120 °C for 1 h, cap the bottle, and cool to room temperature in a desiccator. Momentarily release the cap to equalize the pressure and weigh the capped bottle to the nearest 0.1 mg. Repeat the drying and weighing until there is no further weight loss. Transfer the test sample to a dry 300-mL Erlenmeyer flask and reweigh the capped bottle to the nearest 0.1 mg. The difference is the mass of the test sample.

11.2 Add 50.0 mL of the ferrous ammonium sulfate solution, plus an additional 10.0 mL for each 0.1 g of MnO_2 present, to the flask. Close the flask with a stopper equipped with inlet and outlet tubes. Pass carbon dioxide through the flask.

11.3 Heat the flask moderately and shake intermittently until the ore is decomposed.

11.4 Cool the contents of the flask while maintaining the flow of carbon dioxide.

11.5 Unstopper the flask, add 2 mL of sodium diphenylamine sulfonate indicator solution, and 10 mL of H_3PO_4 . Dilute to 150 mL with cold water (from which the air was removed by boiling) and titrate the excess ferrous ammonium sulfate with standard $\text{K}_2\text{Cr}_2\text{O}_7$ solution to a permanent purple end point.

11.6 The correlation between the solutions of ferrous ammonium sulfate and potassium dichromate is established under test conditions. For this purpose, pour into a flask the same amount of ferrous ammonium sulfate solution used to dissolve the ore, and proceed as directed in 11.3.

12. Calculation

12.1 Calculate the manganese dioxide content as follows:

$$\text{Manganese dioxide, \%} = \frac{[(A - B) \times C \times 4.3465]}{D} \quad (1)$$

where:

A = millilitres of standard $\text{K}_2\text{Cr}_2\text{O}_7$ solution used to establish the correlation in 11.6,

B = millilitres of standard $\text{K}_2\text{Cr}_2\text{O}_7$ solution required to titrate the excess of ferrous ammonium sulfate in the sample solution,

C = the normality of standard $\text{K}_2\text{Cr}_2\text{O}_7$ solution, and

D = grams of test sample used.

NOTE 1—As used in this test method (except as related to the term *relative standard deviation*), “percent” or “%” refers to mass fraction (wt/wt) of the form g/100g.

13. Precision and Bias

13.1 *Precision*—Table 1 indicates the precision of the test method between laboratories.

13.2 *Bias*—No information on the bias of this test method is known. Accepted reference materials may have not been included in the materials used in the interlaboratory study. Users of the method are encouraged to employ accepted reference materials, if available, and to judge the bias of the method from the difference between the accepted value for the manganese dioxide content and the mean value from interlaboratory testing of the reference material.

TEST METHOD B—SODIUM OXALATE METHOD

14. Summary of Test Method

14.1 The test sample is dissolved in sulfuric acid in the presence of sodium oxalate. The manganese dioxide reacts with an equivalent amount of oxalate. The excess sodium oxalate is titrated with a standard solution of potassium permanganate.

15. Reagents and Materials

15.1 *1,10 Phenanthroline Indicator Solution* (0.025 M)

15.1.1 Dissolve 1.485 g of 1,10-phenanthroline monohydrate and 0.695 g of ferrous sulfate ($\text{Fe SO}_4 \cdot 7 \text{H}_2\text{O}$) in 50 mL of water. Dilute to 100 mL.

15.2 *Potassium Permanganate, Standard Solution* (0.1 N)

15.2.1 *Preparation*—Dissolve 3.2 g of potassium permanganate (KMnO_4) in 1 L of water. Let stand in the dark for two weeks. Filter, without washing, through a Gooch crucible or a

TABLE 1 Statistical Information

Sample	Method	Average Manganese Dioxide Content, %	Relative Standard Deviation, %	Number of Determinations	Number of Participating Laboratories
1	A	40.87	0.54	35	7
2	A	70.23	0.47	29	6
1	B	40.90	0.46	29	8
2	B	70.32	0.38	29	8

fine porosity fritted-glass crucible. Avoid contact with rubber or other organic material. Store in a dark-colored glass-stoppered bottle.

15.2.2 Standardization—Dry a portion of the primary standard sodium oxalate at 105 °C. Transfer 0.3000 g of the sodium oxalate to a 600-mL beaker. Add 250 mL of H₂SO₄(5 + 95), previously boiled for 10 to 15 min and then cooled to 27 °C ± 3 °C, and stir until the oxalate has dissolved. Add (39 to 40) mL (Note 2) of the KMnO₄ solution, at a rate of (25 to 35) mL/min, while stirring slowly. Let stand until the pink color disappears (about 45 s) (Note 3). Heat to 55 °C to 60 °C and complete the titration by adding KMnO₄ solution until a faint pink color persists for 30 s. Add the last (0.5 to 1) mL dropwise, allowing each drop to become decolorized before adding the next drop. To determine the blank: titrate 250 mL of H₂SO₄(5 + 95), treated as above, with KMnO₄ solution to a faint pink color. The blank correction is usually equivalent to 0.03 mL to 0.05 mL.

NOTE 2—A 0.3000-g portion of sodium oxalate requires 44.77 mL of KMnO₄ solution (0.1 N) (10.2).

NOTE 3—If the KMnO₄ solution is too strong, the pink color will fade at this point; begin again, adding a few millilitres less of the KMnO₄ solution.

15.3 Sodium Oxalate—Dry the reagent for 2 h at 105 °C prior to use.

16. Procedure

16.1 Transfer approximately 0.4 g of test sample to a small, dry weighing bottle and place in a drying oven. Dry at 120 °C for 1 h, cap the bottle, and cool to room temperature in a desiccator. Momentarily release the cap to equalize pressure and weigh the capped bottle to the nearest 0.1 mg. Repeat the drying and weighing until there is no further weight loss. Transfer the test sample to a dry 300-mL Erlenmeyer flask and reweigh the capped bottle to the nearest 0.1 mg. The difference is the mass of the sample.

16.2 Add 100 mL of H₂SO₄ (1 + 9) and 0.8000 g of sodium oxalate to the flask.

16.3 Cover the flask with a small cover glass and heat on a steam bath to decompose the ore. Swirl the flask occasionally and continue heating until all dark colored particles have disappeared.

16.4 When decomposition is complete, rinse the contents of the flask into a 600-mL beaker and adjust the volume to about 200 mL with hot water. Add 2 drops to 3 drops of the 1,10 phenanthroline indicator solution and titrate the hot solution (60 °C to 70 °C) with the standard KMnO₄ solution (Note 4). At the end point, the color will change from pink to green.

NOTE 4—The titration may be performed without using the indicator by observing a pink end point due to excess potassium permanganate.

16.5 The correlation between the sodium oxalate and the standard permanganate solution is established under test conditions. For this purpose, transfer 0.8000 g of sodium oxalate and 100 mL of H₂SO₄ (1 + 9) to a 300-mL Erlenmeyer flask. Proceed as directed in 16.3.

17. Calculation

17.1 Calculate the manganese dioxide content as follows:

$$\text{Manganese dioxide, \%} = \frac{[(A - B) \times C \times 4.3465]}{D} \quad (2)$$

where:

- A = millilitres of standard KMnO₄ solution used to titrate 0.8000 g of sodium oxalate,
- B = millilitres of standard KMnO₄ solution used to titrate the excess of sodium oxalate in the sample solution,
- C = normality of the standard KMnO₄ solution, and
- D = grams of sample used.

18. Precision and Bias

18.1 Precision—Table 1 indicates the precision of the test method between laboratories.

18.2 Bias—No information on the bias of this test method is known. Accepted reference materials may have not been included in the materials used in the interlaboratory study. Users of the method are encouraged to employ accepted reference materials, if available, and to judge the bias of the method from the difference between the accepted value for the manganese dioxide content of the reference material and the mean value from interlaboratory testing of the reference material.

19. Keywords

19.1 manganese dioxide content; manganese ores; redox titrimetry

ASTM International takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.

This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, at the address shown below.

This standard is copyrighted by ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, United States. Individual reprints (single or multiple copies) of this standard may be obtained by contacting ASTM at the above address or at 610-832-9585 (phone), 610-832-9555 (fax), or service@astm.org (e-mail); or through the ASTM website (www.astm.org). Permission rights to photocopy the standard may also be secured from the ASTM website (www.astm.org/COPYRIGHT/).