



## Standard Specification for Shoulder Prostheses<sup>1</sup>

This standard is issued under the fixed designation F1378; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon ( $\epsilon$ ) indicates an editorial change since the last revision or reapproval.

### 1. Scope

1.1 This specification covers shoulder prostheses for total or hemiarthroplasty used to provide functioning articulation by employing glenoid and humeral components.

1.2 Devices for custom applications are not covered by this specification. Modular prostheses are included in this specification.

1.3 The values stated in SI units are to be regarded as standard. No other units of measurement are included in this standard.

### 2. Referenced Documents

#### 2.1 ASTM Standards:<sup>2</sup>

**F75** Specification for Cobalt-28 Chromium-6 Molybdenum Alloy Castings and Casting Alloy for Surgical Implants (UNS R30075)

**F86** Practice for Surface Preparation and Marking of Metallic Surgical Implants

**F90** Specification for Wrought Cobalt-20Chromium-15Tungsten-10Nickel Alloy for Surgical Implant Applications (UNS R30605)

**F136** Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401)

**F138** Specification for Wrought 18Chromium-14Nickel-2.5Molybdenum Stainless Steel Bar and Wire for Surgical Implants (UNS S31673)

**F562** Specification for Wrought 35Cobalt-35Nickel-20Chromium-10Molybdenum Alloy for Surgical Implant Applications (UNS R30035)

**F563** Specification for Wrought Cobalt-20Nickel-20Chromium-3.5Molybdenum-3.5Tungsten-5Iron Alloy for Surgical Implant Applications (UNS R30563)<sup>3</sup>

**F603** Specification for High-Purity Dense Aluminum Oxide for Medical Application<sup>3</sup>

**F648** Specification for Ultra-High-Molecular-Weight Polyethylene Powder and Fabricated Form for Surgical Implants

**F745** Specification for 18Chromium-12.5Nickel-2.5Molybdenum Stainless Steel for Cast and Solution-Annealed Surgical Implant Applications

**F746** Test Method for Pitting or Crevice Corrosion of Metallic Surgical Implant Materials

**F748** Practice for Selecting Generic Biological Test Methods for Materials and Devices

**F799** Specification for Cobalt-28Chromium-6Molybdenum Alloy Forgings for Surgical Implants (UNS R31537, R31538, R31539)

**F981** Practice for Assessment of Compatibility of Biomaterials for Surgical Implants with Respect to Effect of Materials on Muscle and Bone

**F983** Practice for Permanent Marking of Orthopaedic Implant Components

**F1044** Test Method for Shear Testing of Calcium Phosphate Coatings and Metallic Coatings

**F1108** Specification for Titanium-6Aluminum-4Vanadium Alloy Castings for Surgical Implants (UNS R56406)

**F1147** Test Method for Tension Testing of Calcium Phosphate and Metallic Coatings

**F1537** Specification for Wrought Cobalt-28Chromium-6Molybdenum Alloys for Surgical Implants (UNS R31537, UNS R31538, and UNS R31539)

**F1829** Test Method for Static Evaluation of Glenoid Locking Mechanism in Shear

**F2028** Test Methods for Dynamic Evaluation of Glenoid Loosening or Disassociation

#### 2.2 ANSI Standard:<sup>4</sup>

ASME B46.1–1995

### 3. Terminology

#### 3.1 Definitions of Terms Specific to This Standard:

3.1.1 *collar*—flange at the junction of the neck and stem.

<sup>1</sup> This specification is under the jurisdiction of ASTM Committee F04 on Medical and Surgical Materials and Devices and is the direct responsibility of Subcommittee F04.22 on Arthroplasty.

Current edition approved Jan. 1, 2012. Published January 2012. Originally approved in 1992. Last previous edition approved in 2010 as F1378 – 05 (2010). DOI: 10.1520/F1378-12.

<sup>2</sup> For referenced ASTM standards, visit the ASTM website, [www.astm.org](http://www.astm.org), or contact ASTM Customer Service at [service@astm.org](mailto:service@astm.org). For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

<sup>3</sup> Withdrawn. The last approved version of this historical standard is referenced on [www.astm.org](http://www.astm.org).

<sup>4</sup> Available from American National Standards Institute (ANSI), 25 W. 43rd St., 4th Floor, New York, NY 10036.

3.1.2 *glenoid component*—the prosthetic portion that replaces, in part or in total, the glenoid fossa of the scapula and articulates with the natural humeral head or a prosthetic replacement.

3.1.3 *head*—bearing member for articulation with the glenoid.

3.1.4 *humeral component*—the prosthetic portion that replaces, in part or in toto, the proximal humerus or humeral head and articulates with the natural glenoid fossa or a prosthetic replacement.

3.1.5 *keel, (or pegs)*—single or multiple projections that provide resistance to translation or rotation of the glenoid component, or both, by mating with cavities created in the glenoid fossa.

3.1.6 *neck*—segment connecting the head and the stem.

3.1.7 *reverse design shoulder implants*—implants that have a ball-shaped glenoid component and a concave humeral design.

3.1.8 *stem*—segment intended for insertion within the humeral medullary canal.

## 4. Classification

4.1 *Constrained*—A constrained joint prosthesis is used for joint replacement and resists dislocation of the prosthesis in more than one anatomical plane and consists of either a single, flexible, across-the-joint component or more than one component linked together or affixed.

4.2 *Partially Constrained*—A semi-constrained joint prosthesis is used for partial or total joint replacement and limits translation and rotation of the prosthesis in one or more planes via the geometry of its articulating surfaces. It has no across-the-joint linkages.

4.3 *Unconstrained*—An unconstrained joint prosthesis is used for partial or total joint replacement and restricts minimally prosthesis movement in one or more planes. Its components have no across-the-joint linkage.

## 5. Materials and Manufacture

5.1 The choice of materials is understood to be a necessary but not sufficient assurance of function of the device made from them. All devices conforming to this specification shall be fabricated from materials, with adequate mechanical strength and durability, corrosion resistance, and biocompatibility.

5.1.1 *Mechanical Strength*—Various components of shoulder prostheses have been successfully fabricated from the following materials. However, not all of these materials may possess sufficient mechanical strength for critical highly-stressed components. See Specifications **F75**, **F90**, **F136**, **F138**, **F562**, **F563** (nonbearing use only), **F603**, **F648**, **F745**, **F799**, **F1108**, and **F1537**.

5.1.2 *Corrosion Resistance*—Materials with limited or no history of successful use for orthopedic implant application shall be determined to exhibit corrosion resistance equal to or better than one of the materials listed in **5.1.1** when tested in accordance with Test Method **F746**.

5.1.3 *Biocompatibility*—Materials with limited or no history of successful use for orthopedic implant application must be determined to exhibit acceptable biological response equal

to or better than one of the materials listed in **5.1.1** when tested in accordance with Practices **F748** and **F981**.

## 6. Performance Requirements

6.1 *Wear of Alternative Materials*—It is important to understand the wear performance for articulating surfaces. Any new or different material should not exceed the wear rates of the following material couple when tested under physiological conditions. The current standard wear couple is CoCrMo alloy (Specification **F75**) against ultra high molecular weight polyethylene (Specification **F648**), both having prosthetic quality surface finishes in accordance with **8.2**.

NOTE 1—In situations where the pin-on-flat test may not be considered appropriate, other test methods may be considered.

6.2 *Range of Motion of Shoulder Prosthesis Prior to Implantation*—Flexion shall be equal to or greater than 90°. Abduction shall be equal to or greater than 90°. Internal rotation shall be equal to or greater than 90°. External rotation shall be equal to or greater than 45°. Extension shall be equal to or greater than 45°.

6.3 Porous metal coatings shall be tested according to Test Method **F1044** (shear strength) and Test Method **F1147** (tensile strength).

### 6.4 Guidelines for In-Vitro Laboratory Testing:

6.4.1 Implant testing should reflect current clinical failures and potential failure modes particular to the implant. These tests may be directed towards subluxation, glenoid loosening, insert dissociation from a metal backing, and humeral head dissociation. To facilitate such testing, several references on shoulder forces have been compiled.<sup>5, 6, 7</sup> Based upon the work by Anglin et al<sup>5</sup> and Poppen et al,<sup>6</sup> the normal shoulder joint reaction forces are on the order of 1 to 2 times body weight with the directions of loading being given in Figure 3 of the study by Anglin et al.<sup>5</sup> In the design of shoulder implants, this background information of the forces and their directions may be helpful in determining worst-case shoulder joint forces. However, these joint reaction forces are based upon normal subjects. In order to generate pass/fail criterion (that is, forces, angles, and number of cycles) for a particular shoulder prosthesis, one should take into consideration the anticipated patient population, the worst-case physiological loads and angles, an appropriate safety factor, and the potential for unsupported surfaces.

6.4.2 All modular implants should be tested in accordance with Test Method **F1829**.

6.4.3 All prosthetic glenoid components shall be capable of withstanding sustained static and dynamic physiological forces of up to 1 times body weight (per **6.4.1**) without compromise of their function for the intended use and environment. All implants should be tested for loosening for a clinically relevant

<sup>5</sup> Anglin, C., Wyss, U. P., Pichora, D. R., "Glenohumeral Contact Forces," *Proceedings of the Institution of Mechanical Engineers. Part H—Journal of Engineering in Medicine*, 214 (6), 2000, pp. 637–644.

<sup>6</sup> Poppen, N. K., Walker, P.S., "Forces at the Glenohumeral Joint in Abduction," *Clinical Orthopaedics & Related Research*, 135, Sept. 1978, pp. 165–170.

<sup>7</sup> Anglin, C., Wyss, U. P., Pichora, D. R., "Mechanical Testing of Shoulder Prostheses and Recommendations for Glenoid Design," *J. of Shoulder and Elbow Surgery*, 9, 2000, pp. 323–331.

number of cycles. It has been suggested by Anglin et al<sup>7</sup> that 100 000 cycles is a suitable number of cycles. A larger number of cycles may be required for worst-case situations. One method for testing cemented glenoid components is Test Method F2028.

## 7. Dimensions

7.1 Dimensions of shoulder joint replacement components shall be as designated in Figs. 1-3.

## 8. Finish and Product Marking

8.1 Items conforming to this specification shall be finished and marked in accordance with Practice F86, where applicable.

### 8.2 Articulating Surface Finishes:

8.2.1 *Metallic Bearing Surface*—The main bearing surface shall have a surface finish no rougher than 0.10  $\mu\text{m}$  roughness average,  $R_a$ , with a cutoff length of 0.25 mm, when measured according to the principles given in ASME B46.1–1995.

8.2.2 *Polymeric Bearing Surface (if used)*—The main bearing surface shall have a surface finish no rougher than 2  $\mu\text{m}$  roughness,  $R_a$ , with a cut-off length of 0.8 mm, when measured according to the principles given in ASME B46.1–1995.

8.3 In accordance with Practices F86 and F983, items conforming to this specification shall be marked as follows in order of priority where space permits: manufacturer, material, lot number, catalog number, and size. Additional information may include a designation for left or right and front.

8.3.1 Optional glenoid marking may specify orientation (top, if applicable; right or left, if applicable).

8.4 If one of the components is not radiographically opaque, it is strongly encouraged that it shall contain a marker wire or other means of radiographic detection located at the manufacturer's discretion.


## 9. Labeling

9.1 The dimensions shown in Figs. 1-3 shall be included in the product labeling.

9.2 The material(s) used for the implant shall be specified on the package labels and inserts.


## 10. Keywords

10.1 arthroplasty; glenoid; humeral; prostheses; hemi-shoulder replacement; total shoulder replacement




NOTE—A modular connection may be included in this device.

FIG. 1 Humeral Collarless Design



NOTE—A modular connection may be included in this device.

**FIG. 2 Humeral Collared Design**



NOTE 1—If the glenoid component is not symmetric about the transverse plane, a minimum and maximum component width shall be specified.

NOTE 2—A modular connection may be included in this device.

**FIG. 3 Glenoid Components**

## APPENDIX

### (Nonmandatory Information)

#### X1. RATIONALE

X1.1 The objectives of this specification are the provision of guidelines for the manufacture and use of the components for total shoulder replacement. Total shoulder replacement parts are intended for use in a patient who is skeletally mature, under conditions of imposed dynamic loads, in a corrosive environment, and virtually continuous motion at the bearing surfaces. Laboratory tests to accurately simulate imposed loads, aggressive electrolytes, and complex constituents of

body fluids have not been usefully accelerated at the present time for a complete joint evaluation. Long-term projections of satisfactory performance over many decades can be suggested but not accurately predicted using available screening procedures. This specification identifies those factors considered to be important to assure a satisfactory useful prosthetic life. It is here recognized that failure of an arthroplasty can occur, even while the components are intact. This is true owing to the

composite nature of the arthroplasty procedure, which includes the implant, cement if any, and the physiological environment.

X1.2 Under applicable documents and materials, the list reflects the current state of the art. It is recognized that should materials not now included appear and be proved acceptable, they shall be inserted in the process of revision.

X1.2.1 *Performance Considerations*—Component performance should be considered with regard to patient anatomy. It is well recognized that physical stresses resulting from events or activities out of the ordinary range, as in accidents or especially vigorous sports, predictably exceed allowable stress levels in any component design. It is also recognized here that other forms of arthroplasty failure are known to occur, related primarily to patient factors such as osteoporosis, Paget's disease, misuse and disuse, and others.

X1.2.1.1 No device-specific wear test is specified in this specification. It is felt that at this time wear is not a major issue in existing or potential implant designs, that presently there are

no techniques available to do device-specific wear tests and that this consideration is already partly covered in 6.1.

X1.2.1.2 The range of motion parameters are specified as minimum values so that the implant itself does not restrict patient shoulder motion and thus allow potential dislocation or loosening of the implant. Initially the minimum range of motion parameters were larger but they were decreased to those given in 6.2.

X1.2.2 *Dimensions*—The method of dimensional measurement shall be sought to conform with industry practice and, whenever possible, on an international basis.

X1.2.3 *Finish and Markings*—Dimensions and tolerances are as described by ANSI documents for sphericity, concentricity, and surface finish.

X1.2.3.1 The manufacturer's trademark shall appear legibly on each of the components. It is desirable to have complete information, where space is available to do so, including size, orientation if any, and catalog number with date.

X1.3 This specification was revised in 2004 to include reverse design shoulder implants.

*ASTM International takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.*

*This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, at the address shown below.*

*This standard is copyrighted by ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, United States. Individual reprints (single or multiple copies) of this standard may be obtained by contacting ASTM at the above address or at 610-832-9585 (phone), 610-832-9555 (fax), or service@astm.org (e-mail); or through the ASTM website (www.astm.org). Permission rights to photocopy the standard may also be secured from the ASTM website (www.astm.org/COPYRIGHT/).*